

Hyperbaric - Ozone - Therapy

BIOZON - H. O. T.



Ozon & Colon Medizintechnik

D - 72 805 Lichtenstein

Lilienweg 4

☎ +49-7129-936090 📠 +49-7129-936091

www.biozon.eu • mail@ozon.nu

Hyperbaric - Ozone - Therapy

BIOZON - H. O. T.

with an integrated controller for a blood treatment
offers significant advantages

- ◆ **EASY TO USE:** Simple and safe to handle. The Sequence of a Large blood treatment is easily operated by two press buttons, therefore the therapist can focus on the patient.
- ◆ **Short treatment:** low effort. The Intake - Ozonation and Re-infusion of 150-170ml of blood takes 8-12 minutes.
- ◆ **Further Features:** Syringe filling - Water ozonation - Bag fumigation

<h2>INDICATION</h2>

❖ **PERIPHERAL AND CEREBRAL** circulatory

❖ **RHEUMATIC FORMS** CIRCLE

❖ **Liver damage**

- Metabolic disorders
- Hepatosis
- Coronary insufficiency
- condition after heart attack
- Myocardial damage
- Vegetative dystonia
- Application in the bladder, vagina and urethra
(All infections of these organs)
- Intramuscular, Intra-arterial,
Subcutaneous, intra-articular -injections